



INFRARED

Because of their safe and natural operation, infrared lamps, or heating lamps, are used in physiotherapy, osteopathy, and massage centres to aid post-exercise recovery, stimulate blood flow, reduce pain and eliminate toxins.

Our heating lamp comes mounted on a flexible gooseneck (86 cm long, 34 cm of which are rigid tube). Equipped with a handle it is easily positioned at the required distance. The large opening of the lamp head is covered with a protective grill to avoid the risk of accidental contact with the heat source.

The model comes with a timer switch that enables the length of treatment to be set from zero to 60 minutes with an audible signal upon stopping.



On heavy trolley base (1)



Timer switch option

INFRARED

| | |
|---|------------------------------------|
| Watts | 250 W |
| Wavelength | 1200 nm |
| Bulb lifetime IR (E27) | 5000 hours minimum |
| Mounting accessories | 1 heavy, 3, 4, 5, 6, 7, 10, 11, 12 |
| Cable length | 4 meters |
| Insulation | Class I |
| Available colours | ○ |
| Lamp head size | 260 mm |
| Lamp base size | 69 x 69 x 115 mm |
| Interaxial spacing (4 holes on lamp's base) | 54,5 mm |

CODES

| | On / off switch | Contactless switch | Timer switch |
|---------------------------|-----------------|--------------------|--------------|
| Flexible gooseneck 860 mm | IRT250860 | IRT250860-ISC | IRT250860MI |

